

# BLOCK



### Bone biopsy needle procedure

- 1 - Penetrate the cortical surface
- 2 - Withdraw the bone introducer mandrel
- 3 - Rotate the cutting cannula inside the bone tissue
- 4 - Insert anti luxation B-lock system

### CODES TO ORDER

CODE	DIAMETER (G)	LENGTH (CM)	PACKAGE
VWJT0810	8	10	10
VWJT0815	8	15	20
VWJT0910	9	10	10
VWJT0915	9	15	20
VWJT1110	11	10	10
VWJT1115	11	15	20
VWJT1307	13	07	10
VWJT1310	13	10	10

On request Vigeo is available to provide needles with different diameter and length