

Percutaneous Ethanol Injection Needle (PEIT)



- Larger volume of treated tissue
- Lower ethanol flow into the peritoneal cavity
- No defelction of the needle from the path
- Suitable for small and deep nodules
- Centimeter depth markings assist in needle placement
- Indicator shows flow direction
- Colour coded for quick gauge size identification
- Luer extension tube for a flexible connection to the syringe



CT

ALCOJECT Qty/Box 10pcs		
CODE	SIZE	DESCRIPTION
● VVHA32115	21Gx15cm	Open Tip

On request Vigeo is available to provide needles with different diameter and length