


Periurethral



Transurethral



DURASPHERE - Indicated for the treatment of adults with stress urinary incontinence due to intrinsic sphincteric deficiency (SUI) or to treat vesicoureteral reflux (VUR).

Cost Effective - Significantly lower cost permits more patients to be treated.

Permanent Material - Pyrolytic carbon-coated beads cannot be degraded in the body.

Non-migratory - Injectable beads are larger than the 80-micron threshold for migration.

Biocompatible - Pyrolytic carbon has been used for more than 30 years in mechanical heart valves.

Designed for Surgeons - Delivery system uses thumb pressure on plunger to provide tactile feedback from tissue.


Surgical Choices - Can be injected using Submucosal (modified periurethral) or Transurethral technique.

Convenient - Available in 1 ml and 3 ml syringe sizes.

Non-immunogenic - No skin test required.

DURASPHERE GR - Intended to treat problems associated with Gastroesophageal Reflux Disease (GERD).

DURASPHERE FI - Intended to treat problems associated with Fecal Incontinence in adult men and women.



DURASPHERE
Injectable Bulking Agent

CODE	SYRINGE
030966	1 ml
030967	3 ml (periurethral)

CODE	NEEDLE	DIAMETER (G)
101082	pencil tip (38 cm)	20G
101074	angled tip (3,8 cm)	20G



Axial CT Scan