

RELOCK SOFT WOMEN'S HEALTH


CE
0476

Soft Repositionable Breast Localization Device


PATENT PENDING



- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Retractable hook with unique flexibility of its Polyamid guiding thread
- Indicator on wing shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- 9mm Loop hook



1. Half exposure - 2. Full exposure - 3. Stiffener cannula



Mammography

RELOCK SOFT Qty/Box 10pcs

CODE	SIZE
● VVRPS1905	19Gx5cm
● VVRPS1907	19Gx7cm
● VVRPS1910	19Gx10cm

On request Vigeo is available to provide needles with different diameter and length