

RELOCK PREMIUM DUO WOMEN'S HEALTH

CE
0476

Repositionable One-Handed Breast Localization Device


PATENT PENDING



- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Removable handle for mammography procedure
- Retractable hook with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimeter depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- Double 5mm Loop hook



Half exposure



Full exposure



Removable handle

RELOCK PREMIUM DUO Qty/Box 10 pcs

CODE	SIZE
● VVNSPX1805	18Gx5cm
● VVNSPX1807	18Gx7cm
● VVNSPX1810	18Gx10cm

On request Vigeo is available to provide needles with different diameter and length